



Agilent Technologies

Sample ID: PC3817K AS144323 4-Fluoro-2-methoxyphenol GF Method Name: Apollo_PC_Data_Collection_ATR
Sample Scans: 64 User: Bruce
Background Scans: 64 Date/Time: 15/02/2017 13:07:15
Resolution: 4 cm⁻¹ Range: 4,000.00 - 650.00
System Status: Good Apodization: Triangular
File Location: C:\Program Files\Agilent\MicroLab PC\Results\Apollo\PC Coded compounds\PC3817K AS144323
4-Fluoro-2-methoxyphenol GFA115119_2017-02-15T13-09-42.a2r

